



Colomate

Colomate Mini Stoma Caps

Our **Colomate Mini Stoma Caps** do not have an integral filter like many traditional 'stoma caps'. Instead, they have a **super absorbent pad** designed to absorb moisture emitted from a defunctioned colostomy (mucous fistula) or to absorb drops of urine from a Mitrofanoff continent urinary reservoir if needed in between catheterisations. For support and more information on Mitrofanoffs please visit www.mitrofanoffsupport.co.uk.

Other uses include healing fistulae or sinus drainage which is often preferable to moist wound dressings which require frequent changing.

The benefits to you...

- Extra thin hydrocolloid skin barrier with elastomers – enables greater flexibility, providing a softer more comfortable system.
- No filter.
- Super absorbent pad.
- Beige non-woven both sides – prevents skin irritation by absorbing perspiration.

This product range is latex free.

Available on NHS prescription.

Colomate Mini Stoma Caps Product Codes – Pack of 30

Stoma Size	Diameter	Product Code
Pre-cut 20mm	6cm	CCS120
Pre-cut 25mm	7cm	CCS125
Pre-cut 30mm	8.5cm	CCS130



For more information please call Freephone 0800 652 0424

