



Combimate

## Combimate 2-Piece Ostomy System

The **Combimate 2-Piece Ostomy System** encompasses all the special features and benefits of the 1-piece pouches with the 2-piece user in mind.

### The benefits to you...

#### Combimate 2-Piece Pouches

- A high capacity 32mm bar filter which helps to prevent ballooning and pancaking, contributing to improved pouch wear time.
- Additional layer of film, the 'raincoat', separates the filter from the stoma providing optimum performance.
- Filter covers are provided for additional protection and security, as directed by your Stoma Care Nurse.
- Body side covering is made of a soft, non-woven material that helps to prevent skin soreness by absorbing moisture from the skin.
- Available in beige for discretion or in a 'needle pin' split cover to assist with pouch application and observation of the stoma.
- **All Combimate 2-piece pouches will securely fit our range of Combimate 2-piece flanges.** You can also **Mix 'n' Match** with ConvaTec's Natura\* range of flanges.
- Available in 2-piece closed, drainable and urostomy pouches in a choice of flange sizes.

#### Adjustable Ostomy Belts

- Provide extra comfort and peace of mind when used with our **Combimate 2-Piece Ostomy System** and most other ostomy appliances. Ensuring even distribution of pressure and support around the peristomal area.
- Available in a **25mm** or **50mm** wide belt with 2-point belt lug connection to fit all sizes or in a 50mm wide belt with 4-point belt lug connection in various waist sizes.
- Easy to attach.
- Machine washable.

\*ConvaTec and Natura are registered trademarks of E.R. Squibb and Sons.

**This product range is latex free.**

**Available on NHS prescription.**

For more information please call Freephone 0800 652 0424



## Belt Support Plates

- Provide extra comfort and peace of mind when used with our **Adjustable Ostomy Belts that have a unique 4-point belt attachment feature.**
- Provides even distribution of pressure and support around the peristomal area.
- Available in 3 sizes.
- Easy to attach.

## Combimate 2-Piece Flanges

- The Combimate 2-piece flanges are designed with softness, flexibility and comfort in mind. This has been achieved by laminating the hydrocolloid skin barrier with a flexible thin layer of polymeric foam.
- To provide comfort and protection for the stoma, the flanges have no plastic film which can potentially damage the stoma.
- For mild and severe stomal retraction, we offer a range of **shallow** and **deep** convex flanges\*\*.
- **All Combimate flanges will securely fit our range of Combimate 2-piece pouches.**  
You can also **Mix 'n' Match** with ConvaTec's Natura\* range of pouches.
- Available in corresponding flange gasket sizes: **38mm, 45mm, 57mm** and **70mm**. Also available in a wide selection of pre-cut sizes which can easily be stretched to the exact stoma size.

\*ConvaTec and Natura are registered trademarks of E.R. Squibb and Sons.

\*\*All convex products should be used upon the advice and guidance of a Stoma Care Nurse.

**This product range is latex free.**

**Available on NHS prescription.**

# Combimate 2-Piece Closed Pouches Product Codes

## Large Closed Pouch – Pack of 30

Flange Size	Transparent	Beige
38mm	-	CLB238
45mm	CLT245	CLB245
57mm	CLT257	CLB257
70mm	CLT270	CLB270



## Medium Large Anatomical Closed Pouch – Pack of 30

Flange Size	Transparent	Beige
38mm	-	CMLB238
45mm	-	CMLB245
57mm	-	CMLB257
70mm	-	CMLB270



## Medium Closed Pouch – Pack of 30

Flange Size	Transparent	Beige
38mm	-	CMB238
45mm	-	CMB245
57mm	-	CMB257



## Small Closed Pouch – Pack of 30

Flange Size	Transparent	Beige
45mm	-	CSB245



## Stoma Cap Closed Pouch – Pack of 30

Flange Size	Transparent	Beige
45mm	-	CCB245
57mm	-	CCB257



# Combimate 2-Piece Drainable Pouches Product Codes


## Large Drainable Pouch – Pack of 30

Flange Size	Transparent	Beige
38mm	-	ILBV238
45mm	ILTV245	ILBV245
57mm	ILTV257	ILBV257
70mm	ILTV270	ILBV270



## Curved Beige Clips – Pack of 10

Product Code	Product Code
	CLS10



## Medium Drainable Pouch – Pack of 30

Flange Size	Transparent	Beige
38mm	-	IMBV238
45mm	-	IMBV245
57mm	-	IMBV257



## Small Drainable Pouch – Pack of 30

Flange Size	Transparent	Beige
45mm	-	ISBV245
57mm	-	ISBV257



## Combimate 2-Piece Urostomy Pouches Product Codes

### Large Urostomy Pouch – Pack of 20

#### Flange Size

38mm  
45mm  
57mm  
70mm

#### Transparent

ULT238B  
ULT245B  
ULT257B  
ULT270B



### Medium Urostomy Pouch – Pack of 20

#### Flange Size

38mm  
45mm  
57mm

#### Transparent

UMT238B  
UMT245B  
UMT257B

## Urostomy Night Drainage Adaptors Product Code – Pack of 10

One Size

AD10



## Night Drainage Bag Product Code – Pack of 10

One Size

NDB10



## Adjustable Ostomy Belts

### Product Codes – Pack of 1

Belt Width	Waist Size	Product Code
25mm*	One size fits all	SB25
50mm**	One size fits all	SB52
50mm***	Small (60-75cm)	SB50S
50mm***	Medium (75-90cm)	SB50M
50mm***	Large (90-110cm)	SB50L
50mm***	Extra Large (110-130cm)	SB50XL
50mm***	Extra Extra Large (130cm and above)	SB50XXL

\*With 2-point belt lug connection with the option to use a 2 to 1 adaptor to allow the use of the 4-point belt attachment when used with a belt support plate.

\*\*With 2-point belt lug connection.

\*\*\*With 4-point belt lug connection when used with a belt support plate.



25mm – SB25




50mm – SB52



## Belt Support Plates

### Product Codes – Pack of 1

Size	Product Code
60mm	FP60
72mm	FP72
80mm	FP80



For more information please call Freephone 0800 652 0424

## Combimate 2-Piece Flanges Product Codes

### SuperFlex Flange – Pack of 5

Stoma Size	38mm Flange	45mm Flange	57mm Flange	70mm Flange
13mm starter hole	FFB3813	FFB4513	FFB5713	FFB7013
Pre-cut 22mm	FFB3822	FFB4522	-	-
Pre-cut 25mm	FFB3825	FFB4525	-	-
Pre-cut 28mm	FFB3828	FFB4528	-	-
Pre-cut 32mm	-	FFB4532	FFB5732	-
Pre-cut 35mm	-	FFB4535	FFB5735	-
Pre-cut 38mm	-	-	FFB5738	-



### SuperFlex Flange – Pack of 5

Stoma Size	38mm Flange	45mm Flange	57mm Flange	70mm Flange
13mm starter hole	FSS3813	FSS4513	FSS5713	FSS7013
Pre-cut 22mm	FSS3822	FSS4522	-	-
Pre-cut 25mm	FSS3825	FSS4525	-	-
Pre-cut 28mm	FSS3828	FSS4528	-	-
Pre-cut 32mm	-	FSS4532	FSS5732	-
Pre-cut 35mm	-	FSS4535	FSS5735	-
Pre-cut 38mm	-	-	FSS5738	-



### SuperFlex Flange With Shallow Convexity – Pack of 5

Stoma Size	38mm Flange	45mm Flange	57mm Flange	70mm Flange
13mm starter hole	FC3813	FC4513	FC5713	FC7013
Pre-cut 16mm	FC3816	-	-	-
Pre-cut 19mm	FC3819	-	-	-
Pre-cut 22mm	-	FC4522	-	-
Pre-cut 25mm	-	FC4525	-	-
Pre-cut 28mm	-	FC4528	-	-
Pre-cut 32mm	-	-	FC5732	-
Pre-cut 35mm	-	-	FC5735	-
Pre-cut 38mm	-	-	FC5738	-
Pre-cut 41mm	-	-	FC5741	-



### SuperFlex Flange With Deep Convexity – Pack of 5

Stoma Size	45mm Flange	57mm Flange	57mm Flange
13mm starter hole	FXC4513	FXC5713	FXC7013
Pre-cut 19mm	FXC4519	-	-
Pre-cut 22mm	FXC4522	-	-
Pre-cut 25mm	-	FXC5725	-
Pre-cut 28mm	-	FXC5728	-
Pre-cut 32mm	-	FXC5732	FXC7032
Pre-cut 35mm	-	-	FXC7035
Pre-cut 39mm	-	-	FXC7039
Pre-cut 44mm	-	-	FXC7044



**NEW**

## Combimate Soft Convex Flanges Product Codes

### Soft Convex Flange – Pack of 5

**Stoma Size**

13mm start hole  
Pre-cut 16mm  
Pre-cut 19mm  
Pre-cut 22mm  
Pre-cut 25mm  
Pre-cut 30mm  
Pre-cut 35mm  
Pre-cut 41mm

**57mm Flange**

FSC5713-36  
FSC5716  
FSC5719  
FSC5722  
FSC5725  
FSC5730  
FSC5735  
FSC5741

