



Varimate Resist 'Eye Shaped' Washer

The **Varimate 'Eye Shaped' Washer** has the same unique features and benefits as the Varimate Resist Washers. The 'eye shaped' washer helps to fill large irregularities around the stoma.

The benefits to you...

- Designed to help improve pouch wear time, the washer can be moulded to form a thicker seal on either side of the stoma to help prevent unwanted leaks and potential skin breakdown.
- The mouldable washer can be increased in size, and will remain in the desired shape without breaking.
- As the name suggests, the washers are resistant to the effects of urine as well as the aggressive output from faecal stomas.
- The washer remains sticky when rolled and stretched repeatedly, and also when wet.

Available on NHS prescription.

Varimate Resist 'Eye Shaped' Washer Product Code – Pack of 10

Diameter
18-35mm starter hole

Product Code
VMR EYE



For more information please call Freephone 0800 652 0424

