



UltraSmooth Meatal Dilatation Catheter

The **UltraSmooth Meatal Dilatation Catheter** is 40cm long and is a single use hydrophilic coated catheter without drainage eyes, designed specifically for dilatation of urethral strictures. Available in various Charriere sizes to suit all needs.

The benefits to you...

- Durable, consistently stable and highly lubricous, the hydrophilic coating ensures ultra-smooth insertion and removal.
- The durable hydrophilic coating resists 'drying out', unlike other hydrophilic coatings.
- The coating is activated in water within 15 seconds.
- Finger holes allow primary packaging to be opened easily.
- Integrated funnel for grip and no-touch insertion.
- An adhesive pad on the packaging allows it to stick securely to most surfaces for preparation.
- Foil packaging ensures that the product reaches the end user in optimum condition.

Available on NHS prescription.

UltraSmooth Meatal Dilatation Catheter Product Codes

Dilatation Catheter Male Length 40cm – Pack of 25

Charriere Size	Product Code
12Ch	USD4012
14Ch	USD4014
16Ch	USD4016
18Ch	USD4018



For more information please call Freephone 0800 652 0424

